

Antibodies Product Datasheet

Product Name:	E-Tag Antibody
Catalog Number:	C40011
Reactivity:	Species independent
Applications:	WB, ELISA
Dilutions:	WB 1:1000-5000, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic peptide: GAPVPYPDPLEPR conjugated to KLH.
Concentration:	1 mg/ml
Gene Name:	E Tag
Protein Name:	E Tag
Other Name:	E epitope tag antibody GAPVPYPDPLEPR epitope tag antibody GAPVPYPDPLEPR tag antibody
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	Rabbit polyclonal to E tag
Storage and Stability:	-20°C/1 year