

Antibodies Product Datasheet

Product Name:	GAPDH Antibody
Catalog Number:	C40021
Reactivity:	Human, Mouse, Rat, Rb, Ch, Mk, sheep, X
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:5000 IHC: 1:200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of GAPDH
Concentration:	
Gene Name:	GAPDH
Protein Name:	Glyceraldehyde-3-phosphate dehydrogenase
Other Name:	GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH; YM3215
Human Swiss Prot No.:	P04406
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous GAPDH protein.
Storage and Stability:	-20°C/1 year