

## Antibodies Product Datasheet

|                               |                                                                                                               |
|-------------------------------|---------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Plant GST Rabbbit Antibody                                                                                    |
| <b>Catalog Number:</b>        | C40024                                                                                                        |
| <b>Reactivity:</b>            | Nicotianatabacum                                                                                              |
| <b>Applications:</b>          | WB                                                                                                            |
| <b>Dilutions:</b>             | WB: 1:3000-10000                                                                                              |
| <b>Modification:</b>          | Unmodified/Total                                                                                              |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                             |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen. |
| <b>Immunogen:</b>             | Recombinant Protein of Plant GST                                                                              |
| <b>Concentration:</b>         |                                                                                                               |
| <b>Gene Name:</b>             |                                                                                                               |
| <b>Protein Name:</b>          |                                                                                                               |
| <b>Other Name:</b>            | YM3219                                                                                                        |
| <b>Human Swiss Prot No.:</b>  |                                                                                                               |
| <b>Formulation:</b>           | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.                         |
| <b>Specificity:</b>           | The antibody detects endogenous Plant GST proteins.                                                           |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                  |