

Antibodies Product Datasheet

Product Name:	Integrin α IIb Antibody
Catalog Number:	C40027
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Integrin α IIb
Concentration:	1 mg/ml
Gene Name:	ITGA2B GP2B ITGAB
Protein Name:	Integrin α IIb
Other Name:	Integrin alpha-IIb (GPalpha IIb) (GPIIb) (Platelet membrane glycoprotein IIb) (CD antigen CD41) [Cleaved into: Integrin alpha-IIb heavy chain; Integrin alpha-IIb light chain, form 1; Integrin alpha-IIb light chain, form 2]
Human Swiss Prot No.:	P08514
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Integrin α IIb
Storage and Stability:	-20°C/1 year