

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                                                           |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Cdk9 Antibody                                                                                                                                                                                                                             |
| <b>Catalog Number:</b>        | C40036                                                                                                                                                                                                                                    |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                                                                                                                                         |
| <b>Applications:</b>          | WB, ELISA                                                                                                                                                                                                                                 |
| <b>Dilutions:</b>             | WB 1:1000-2000 ELISA 1:5000-20000                                                                                                                                                                                                         |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                                                                                          |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                                                                                         |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                                                                                                 |
| <b>Immunogen:</b>             | Synthesized peptide derived from human Cdk9                                                                                                                                                                                               |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                                                                                   |
| <b>Gene Name:</b>             | CDK9 CDC2L4 TAK                                                                                                                                                                                                                           |
| <b>Protein Name:</b>          | Cdk9                                                                                                                                                                                                                                      |
| <b>Other Name:</b>            | Cyclin-dependent kinase 9 (EC 2.7.11.22) (EC 2.7.11.23) (C-2K) (Cell division cycle 2-like protein kinase 4) (Cell division protein kinase 9) (Serine/threonine-protein kinase PITALRE) (Tat-associated kinase complex catalytic subunit) |
| <b>Human Swiss Prot No.:</b>  | P50750                                                                                                                                                                                                                                    |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                   |
| <b>Specificity:</b>           | This antibody detects endogenous levels of Human Cdk9                                                                                                                                                                                     |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                                                                              |