

Antibodies Product Datasheet

Product Name:	DCT Antibody
Catalog Number:	C40042
Reactivity:	Human, Mouse, (93%)
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 171-244
Concentration:	1 mg/ml
Gene Name:	DCT TYRP2
Protein Name:	DCT
Other Name:	L-dopachrome tautomerase (DCT) (DT) (EC 5.3.3.12) (L-dopachrome Delta-isomerase) (Tyrosinase-related protein 2) (TRP-2) (TRP2)
Human Swiss Prot No.:	P40126
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous DCT
Storage and Stability:	-20°C/1 year