

Antibodies Product Datasheet

Product Name:	CACNB3 Antibody
Catalog Number:	C40044
Reactivity:	Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:500-1000 IHC: 1:200-500
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of CACNB3
Concentration:	
Gene Name:	CACNB3
Protein Name:	Voltage-dependent L-type calcium channel subunit beta-3
Other Name:	CACNB3; CACNLB3; Voltage-dependent L-type calcium channel subunit beta-3; CAB3; Calcium channel voltage-dependent subunit beta 3; YM3398
Human Swiss Prot No.:	P54284
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CACNB3 protein.
Storage and Stability:	-20°C/1 year