

Antibodies Product Datasheet

Product Name:	Cav3.2 Antibody
Catalog Number:	C40049
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC 1:50-100
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Cav3.2
Concentration:	1 mg/ml
Gene Name:	CACNA1H
Protein Name:	Voltage-dependent T-type calcium channel subunit alpha-1H (Low-voltage-activated calcium channel alpha1 3.2 subunit) (Voltage-gated calcium channel subunit alpha Cav3.2)
Other Name:	CACNA1H; Voltage-dependent T-type calcium channel subunit alpha-1H; Low-voltage-activated calcium channel alpha1 3.2 subunit; Voltage-gated calcium channel subunit alpha Cav3.2; YM3537
Human Swiss Prot No.:	O95180
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cav3.2 protein(A208) detects endogenous levels of Cav3.2
Storage and Stability:	-20°C/1 year