

Antibodies Product Datasheet

Product Name:	GABA A Receptor α 4 Antibody
Catalog Number:	C40055
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of GABA A Receptor α 4
Concentration:	1 mg/ml
Gene Name:	GABRA4
Protein Name:	Gamma-aminobutyric acid receptor subunit alpha-4 (GABA(A) receptor subunit alpha-4)
Other Name:	GABRA4; Gamma-aminobutyric acid receptor subunit alpha-4; GABA(A) receptor subunit alpha-4; YM3552
Human Swiss Prot No.:	P48169
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABA A Receptor α 4 protein(A226) detects endogenous levels of GABA A Receptor α 4
Storage and Stability:	-20°C/1 year