

Antibodies Product Datasheet

Product Name:	GABA B Receptor 2 Antibody
Catalog Number:	C40057
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of GABA B Receptor 2
Concentration:	1 mg/ml
Gene Name:	GABBR2
Protein Name:	Gamma-aminobutyric acid type B receptor subunit 2 (GABA-B receptor 2) (GABA-B-R2) (GABA-BR2) (GABABR2) (Gb2) (G-protein coupled receptor 51) (HG20)
Other Name:	GABBR2; GPR51; GPRC3B; Gamma-aminobutyric acid type B receptor subunit 2; GABA-B receptor 2; GABA-B-R2; GABA-BR2; GABABR2; Gb2; G-protein coupled receptor 51; HG20; YM3553
Human Swiss Prot No.:	O75899
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABA B Receptor 2 protein(A228) detects endogenous levels of GABA B Receptor 2
Storage and Stability:	-20°C/1 year