

Antibodies Product Datasheet

Product Name:	KChIP3 Antibody
Catalog Number:	C40060
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of KChIP3
Concentration:	1 mg/ml
Gene Name:	KCNIP3
Protein Name:	Calsenilin (A-type potassium channel modulatory protein 3) (DRE-antagonist modulator) (DREAM) (Kv channel-interacting protein 3) (KChIP3)
Other Name:	KCNIP3; CSEN; DREAM; KCHIP3; Calsenilin; A-type potassium channel modulatory protein 3; DRE-antagonist modulator; DREAM; Kv channel-interacting protein 3; KChIP3; YM3570
Human Swiss Prot No.:	Q9Y2W7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KChIP3 protein(A250) detects endogenous levels of KChIP3
Storage and Stability:	-20°C/1 year