

Antibodies Product Datasheet

Product Name:	Kv10.1 Antibody
Catalog Number:	C40068
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Kv10.1
Concentration:	1 mg/ml
Gene Name:	KCNG3
Protein Name:	Potassium voltage-gated channel subfamily G member 3 (Voltage-gated potassium channel subunit Kv10.1) (Voltage-gated potassium channel subunit Kv6.3)
Other Name:	KCNG3; Potassium voltage-gated channel subfamily G member 3; Voltage-gated potassium channel subunit Kv10.1; Voltage-gated potassium channel subunit Kv6.3; YM3576
Human Swiss Prot No.:	Q8TAE7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Kv10.1 protein(A259) detects endogenous levels of Kv10.1
Storage and Stability:	-20°C/1 year