

Antibodies Product Datasheet

Product Name:	Kv11.3 Antibody
Catalog Number:	C40070
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Kv11.3
Concentration:	1 mg/ml
Gene Name:	KCNH7
Protein Name:	Potassium voltage-gated channel subfamily H member 7 (Ether-a-go-go-related gene potassium channel 3) (ERG-3) (Eag-related protein 3) (Ether-a-go-go-related protein 3) (hERG-3) (Voltage-gated potassium channel subunit Kv11.3)
Other Name:	Potassium voltage-gated channel subfamily H member 7 (Ether-a-go-go-related gene potassium channel 3) (ERG-3) (Eag-related protein 3) (Ether-a-go-go-related protein 3) (hERG-3) (Voltage-gated potassium channel subunit Kv11.3); YM3579
Human Swiss Prot No.:	Q9NS40
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Kv11.3 protein(A263) detects endogenous levels of Kv11.3
Storage and Stability:	-20°C/1 year