

Antibodies Product Datasheet

Product Name:	Flotillin-1 Antibody
Catalog Number:	C40073
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Flotillin-1
Concentration:	1 mg/ml
Gene Name:	FLOT1
Protein Name:	Flotillin-1
Other Name:	Flotillin-1
Human Swiss Prot No.:	O75955
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Flotillin-1
Storage and Stability:	-20°C/1 year