

Antibodies Product Datasheet

Product Name:	HSC 70 Antibody
Catalog Number:	C40079
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:50-300. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human HSC 70.
Concentration:	
Gene Name:	HSPA8
Protein Name:	Heat shock cognate 71 kDa protein
Other Name:	HSPA8; HSC70; HSP73; HSPA10; Heat shock cognate 71 kDa protein; Heat shock 70 kDa protein 8; YM3359
Human Swiss Prot No.:	P11142
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous HSC 70 protein.
Storage and Stability:	-20°C/1 year