

Antibodies Product Datasheet

Product Name:	ANXA5 Antibody
Catalog Number:	C40082
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 60-140
Concentration:	1 mg/ml
Gene Name:	ANXA5 ANX5 ENX2 PP4
Protein Name:	Annexin A5 (Anchored CII) (Annexin V) (Annexin-5) (Calphobindin I) (CBP-I) (Endonexin II) (Lipocortin V) (Placental anticoagulant protein 4) (PP4) (Placental anticoagulant protein I) (PAP-I) (Thromboplastin inhibitor) (Vascular anticoagulant-alpha) (VAC-a
Other Name:	
Human Swiss Prot No.:	P08758
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ANXA5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year