

Antibodies Product Datasheet

Product Name:	CD82 Antibody
Catalog Number:	C40083
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 120-200
Concentration:	1 mg/ml
Gene Name:	CD82 KAI1 SAR2 ST6 TSPAN27
Protein Name:	CD82 antigen (C33 antigen) (IA4) (Inducible membrane protein R2) (Metastasis suppressor Kangai-1) (Suppressor of tumorigenicity 6 protein) (Tetraspanin-27) (Tspan-27) (CD antigen CD82)
Other Name:	
Human Swiss Prot No.:	P27701
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CD82 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year