

Antibodies Product Datasheet

Product Name:	Centaurin- β 1 Antibody
Catalog Number:	C40084
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Centaurin- β 1
Concentration:	1 mg/ml
Gene Name:	ACAP1 CENTB1 KIAA0050
Protein Name:	Centaurin- β 1
Other Name:	Arf-GAP with coiled-coil, ANK repeat and PH domain-containing protein 1 (Centaurin-beta-1) (Cnt-b1)
Human Swiss Prot No.:	Q15027
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Centaurin- β 1
Storage and Stability:	-20°C/1 year