

Antibodies Product Datasheet

Product Name:	CCL17 Antibody
Catalog Number:	C40090
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CCL17
Concentration:	1 mg/ml
Gene Name:	CCL17 SCYA17 TARC
Protein Name:	CCL17
Other Name:	C-C motif chemokine 17 (CC chemokine TARC) (Small-inducible cytokine A17) (Thymus and activation-regulated chemokine)
Human Swiss Prot No.:	Q92583
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CCL17
Storage and Stability:	-20°C/1 year