

Antibodies Product Datasheet

Product Name:	Calponin 1 Antibody
Catalog Number:	C40105
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Calponin 1
Concentration:	1 mg/ml
Gene Name:	CNN1
Protein Name:	Calponin 1
Other Name:	Calponin-1 (Basic calponin) (Calponin H1, smooth muscle)
Human Swiss Prot No.:	P51911
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Calponin 1
Storage and Stability:	-20°C/1 year