

Antibodies Product Datasheet

Product Name:	CD43 Antibody
Catalog Number:	C40106
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CD43
Concentration:	1 mg/ml
Gene Name:	SPN CD43
Protein Name:	CD43
Other Name:	Leukosialin (Galactoglycoprotein) (GALGP) (Leukocyte sialoglycoprotein) (Sialophorin) (CD antigen CD43)
Human Swiss Prot No.:	P16150
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CD43
Storage and Stability:	-20°C/1 year