

Antibodies Product Datasheet

Product Name:	AMPK β 2 Antibody
Catalog Number:	C40118
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AMPK β 2
Concentration:	1 mg/ml
Gene Name:	PRKAB2
Protein Name:	AMPK β 2
Other Name:	5'-AMP-activated protein kinase subunit beta-2 (AMPK subunit beta-2)
Human Swiss Prot No.:	O43741
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human AMPK β 2
Storage and Stability:	-20°C/1 year