

Antibodies Product Datasheet

Product Name:	DUSP6 Antibody
Catalog Number:	C40119
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-3000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of DUSP6
Concentration:	1 mg/ml
Gene Name:	DUSP6
Protein Name:	Dual specificity protein phosphatase 6 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase PYST1) (Mitogen-activated protein kinase phosphatase 3) (MAP kinase phosphatase 3) (MKP-3)
Other Name:	DUSP6; MKP3; PYST1; Dual specificity protein phosphatase 6; Dual specificity protein phosphatase PYST1; Mitogen-activated protein kinase phosphatase 3; MAP kinase phosphatase 3; MKP-3; YM3580
Human Swiss Prot No.:	Q16828
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DUSP6 Rabbit Polyclonal detects endogenous levels of DUSP6
Storage and Stability:	-20°C/1 year