

Antibodies Product Datasheet

Product Name:	ERK 5 Antibody
Catalog Number:	C40139
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:1000-2000 IHC:1:200-500
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of ERK 5
Concentration:	
Gene Name:	MAPK7,ERK5
Protein Name:	Mitogen-activated protein kinase 7
Other Name:	MAPK7; BMK1; ERK5; PRKM7; Mitogen-activated protein kinase 7; MAP kinase 7; MAPK 7; Big MAP kinase 1; BMK-1; Extracellular signal-regulated kinase 5; ERK-5; YM3369
Human Swiss Prot No.:	Q13164
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous ERK 5 protein.
Storage and Stability:	-20°C/1 year