

Antibodies Product Datasheet

Product Name:	AMPK α 2 Antibody
Catalog Number:	C40147
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of AMPK α 2
Concentration:	
Gene Name:	PRKAA2
Protein Name:	5'-AMP-activated protein kinase catalytic subunit alpha-2
Other Name:	PRKAA2; AMPK; AMPK2; 5'-AMP-activated protein kinase catalytic subunit alpha-2; AMPK subunit alpha-2; Acetyl-CoA carboxylase kinase; ACACA kinase; Hydroxymethylglutaryl-CoA reductase kinase; HMGCR kinase; YM3362
Human Swiss Prot No.:	P54646
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous AMPK α 2 protein.
Storage and Stability:	-20°C/1 year