

Antibodies Product Datasheet

Product Name:	AMPK β 2 Antibody
Catalog Number:	C40148
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of AMPK β 2
Concentration:	
Gene Name:	PRKAB2
Protein Name:	5'-AMP-activated protein kinase subunit beta-2
Other Name:	PRKAB2; 5'-AMP-activated protein kinase subunit beta-2; AMPK subunit beta-2; YM3447
Human Swiss Prot No.:	O43741
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous AMPK β 2 protein.
Storage and Stability:	-20°C/1 year