

Antibodies Product Datasheet

Product Name:	ATM Rabbit Antibody
Catalog Number:	C40150
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of ATM of ATM
Concentration:	1 mg/ml
Gene Name:	ATM
Protein Name:	ATM
Other Name:	ATM; YM3684
Human Swiss Prot No.:	Q13315
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATM protein detects endogenous levels of ATM
Storage and Stability:	-20°C/1 year