

Antibodies Product Datasheet

Product Name:	BEST2 Antibody
Catalog Number:	C40151
Reactivity:	Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:500-1000 IHC: 1:200-500
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic Peptide of BEST2
Concentration:	
Gene Name:	BEST2
Protein Name:	Bestrophin-2
Other Name:	Bestrophin-2; Vitelliform macular dystrophy 2-like protein 1; YM3384
Human Swiss Prot No.:	Q8NFU1
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous BEST2 protein.
Storage and Stability:	-20°C/1 year