

Antibodies Product Datasheet

Product Name:	JNK2 Rabbit Antibody
Catalog Number:	C40161
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:500-2000,IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of JNK2
Concentration:	1 mg/ml
Gene Name:	MAPK9 JNK2 PRKM9 SAPK1A
Protein Name:	Mitogen-activated protein kinase 9 (MAP kinase 9) (MAPK 9) (EC 2.7.11.24) (JNK-55) (Stress-activated protein kinase 1a) (SAPK1a) (Stress-activated protein kinase JNK2) (c-Jun N-terminal kinase 2)
Other Name:	Mitogen-activated protein kinase 9 (MAP kinase 9) (MAPK 9) (EC 2.7.11.24) (JNK-55) (Stress-activated protein kinase 1a) (SAPK1a) (Stress-activated protein kinase JNK2) (c-Jun N-terminal kinase 2); YM3721
Human Swiss Prot No.:	P45984
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous JNK2 protein
Storage and Stability:	-20°C/1 year