

Antibodies Product Datasheet

Product Name:	MEK-5 Antibody
Catalog Number:	C40167
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:1000-2000 IHC: 1:200-500
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of MEK-5
Concentration:	
Gene Name:	MAP2K5
Protein Name:	Dual specificity mitogen-activated protein kinase kinase 5
Other Name:	MAP2K5; MEK5; MKK5; PRKMK5; Dual specificity mitogen-activated protein kinase kinase 5; MAP kinase kinase 5; MAPKK 5; MAPK/ERK kinase 5; MEK 5; YM3436
Human Swiss Prot No.:	Q13163
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous MEK-5 protein.
Storage and Stability:	-20°C/1 year