

Antibodies Product Datasheet

Product Name:	CD86 Antibody
Catalog Number:	C40175
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CD86
Concentration:	1 mg/ml
Gene Name:	CD86 CD28LG2
Protein Name:	CD86
Other Name:	T-lymphocyte activation antigen CD86 (Activation B7-2 antigen) (B70) (BU63) (CTLA-4 counter-receptor B7.2) (FUN-1) (CD antigen CD86) P42081
Human Swiss Prot No.:	P42081
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human CD86
Storage and Stability:	-20°C/1 year