

Antibodies Product Datasheet

Product Name:	BMPR-II Antibody
Catalog Number:	C40179
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BMPR-II
Concentration:	1 mg/ml
Gene Name:	BMPR2 PPH1
Protein Name:	BMPR-II
Other Name:	Bone morphogenetic protein receptor type-2 (BMP type-2 receptor) (BMPR-2) (EC 2.7.11.30) (Bone morphogenetic protein receptor type II) (BMP type II receptor) (BMPR-II)
Human Swiss Prot No.:	Q13873
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat, Monkey BMPR-II
Storage and Stability:	-20°C/1 year