

Antibodies Product Datasheet

Product Name:	DLG4 Antibody
Catalog Number:	C40197
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DLG4
Concentration:	1 mg/ml
Gene Name:	DLG4 PSD95
Protein Name:	DLG4
Other Name:	Disks large homolog 4 (Postsynaptic density protein 95) (PSD-95) (Synapse-associated protein 90) (SAP-90) (SAP90)
Human Swiss Prot No.:	P78352
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat DLG4
Storage and Stability:	-20°C/1 year