

Antibodies Product Datasheet

Product Name:	EAA2 Antibody
Catalog Number:	C40198
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EAA2
Concentration:	1 mg/ml
Gene Name:	SLC1A2 EAAT2 GLT1
Protein Name:	EAA2
Other Name:	Excitatory amino acid transporter 2 (Glutamate/aspartate transporter II) (Sodium-dependent glutamate/aspartate transporter 2) (Solute carrier family 1 member 2)
Human Swiss Prot No.:	P43004
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human,Rat,Mouse EAA2
Storage and Stability:	-20°C/1 year