

Antibodies Product Datasheet

Product Name:	EAAT2 Antibody
Catalog Number:	C40200
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of EAAT2
Concentration:	1 mg/ml
Gene Name:	SLC1A2
Protein Name:	Excitatory amino acid transporter 2 (Glutamate/aspartate transporter II) (Sodium-dependent glutamate/aspartate transporter 2) (Solute carrier family 1 member 2)
Other Name:	Excitatory amino acid transporter 2 (Glutamate/aspartate transporter II) (Sodium-dependent glutamate/aspartate transporter 2) (Solute carrier family 1 member 2); YM3549
Human Swiss Prot No.:	P43004
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EAAT2 protein(A223) detects endogenous levels of EAAT2
Storage and Stability:	-20°C/1 year