

Antibodies Product Datasheet

Product Name:	GABA Transporter 1 Antibody
Catalog Number:	C40201
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of GABA Transporter 1
Concentration:	1 mg/ml
Gene Name:	SLC6A1
Protein Name:	Sodium- and chloride-dependent GABA transporter 1 (GAT-1) (Solute carrier family 6 member 1)
Other Name:	SLC6A1; GABATR; GABT1; GAT1; Sodium- and chloride-dependent GABA transporter 1; GAT-1; Solute carrier family 6 member 1; YM3555
Human Swiss Prot No.:	P30531
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GABA Transporter 1 protein(A230) detects endogenous levels of GABA Transporter 1
Storage and Stability:	-20°C/1 year