

Antibodies Product Datasheet

Product Name:	GAP43 Antibody
Catalog Number:	C40202
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:1000-3000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Neuromodulin
Concentration:	
Gene Name:	GAP43
Protein Name:	Neuromodulin
Other Name:	GAP43; Neuromodulin; Axonal membrane protein GAP-43; Growth-associated protein 43; Neural phosphoprotein B-50; pp46; YM3250 P17677
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous GAP43 protein.
Storage and Stability:	-20°C/1 year