

Antibodies Product Datasheet

Product Name:	GSDME n-terminal Antibody
Catalog Number:	C40205
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	IF: 1:50-200 WB: 1:5000-20000 IHC: 1:200-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GSDME. AA range:1-100
Concentration:	1 mg/ml
Gene Name:	GSDME
Protein Name:	Gasdermin- E,DFNA5,Dfna5h,EG14210,2310037D07Rik .4932441K13Rik ,ICERE 1,ICERE-1,Deafness, autosomal dominant 5,Non-syndromic hearing impairment protein 5,GSDME-NT
Other Name:	
Human Swiss Prot No.:	O60443
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	The antibody detects endogenous full length and n-terminal fragment of gsdme protein,
Storage and Stability:	-20°C/1 year