

Antibodies Product Datasheet

Product Name:	IL-3 Antibody
Catalog Number:	C40207
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IL-3
Concentration:	1 mg/ml
Gene Name:	IL3
Protein Name:	IL-3
Other Name:	Interleukin-3 (IL-3) (Hematopoietic growth factor) (Mast cell growth factor) (MCGF) (Multipotential colony-stimulating factor) (P-cell-stimulating factor)
Human Swiss Prot No.:	P08700
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat IL-3
Storage and Stability:	-20°C/1 year