

Antibodies Product Datasheet

Product Name:	Collagen I Antibody
Catalog Number:	C40218
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Collagen I . at AA range: 1200-1217
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	Collagen I
Other Name:	collagen, type I, alpha 2
Human Swiss Prot No.:	P08123/P02452
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Collagen I Polyclonal Antibody detects endogenous levels of Collagen I
Storage and Stability:	-20°C/1 year