

Antibodies Product Datasheet

Product Name:	Adrenocorticotropin(ACTH) Antibody
Catalog Number:	C40259
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Adrenocorticotropin(ACTH)
Concentration:	1 mg/ml
Gene Name:	POMC
Protein Name:	Adrenocorticotropin(ACTH)
Other Name:	Pro-opiomelanocortin (POMC) (Corticotropin-lipotropin) [Cleaved into: NPP; Melanotropin gamma (Gamma-MSH); Potential peptide; Corticotropin (Adrenocorticotropic hormone) (ACTH); Melanotropin alpha (Alpha-MSH); Corticotropin-like intermediary peptide (CLIP
Human Swiss Prot No.:	P01189
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Adrenocorticotropin(ACTH)
Storage and Stability:	-20°C/1 year