

Antibodies Product Datasheet

Product Name:	APOC1 Antibody
Catalog Number:	C40307
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human APOC1
Concentration:	1 mg/ml
Gene Name:	APOC1 APOC1B
Protein Name:	APOC1
Other Name:	Apolipoprotein C-I (Apo-CIB) (ApoC-IB) (Apolipoprotein C1) [Cleaved into: Truncated apolipoprotein C-I (Apo-CIB') (ApoC-IB')]
Human Swiss Prot No.:	P02654
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human APOC1
Storage and Stability:	-20°C/1 year