

Antibodies Product Datasheet

Product Name:	APOC2 Antibody
Catalog Number:	C40308
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human APOC2
Concentration:	1 mg/ml
Gene Name:	APOC2 APC2
Protein Name:	APOC2
Other Name:	Apolipoprotein C-II (Apo-CII) (ApoC-II) (Apolipoprotein C2)
Human Swiss Prot No.:	P02655
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human APOC2
Storage and Stability:	-20°C/1 year