

Antibodies Product Datasheet

Product Name:	BCE1 Antibody
Catalog Number:	C40346
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1: 50-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BCE1
Concentration:	1 mg/ml
Gene Name:	BCE1
Protein Name:	BCE1
Other Name:	
Human Swiss Prot No.:	O60756
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BCE1 at Human
Storage and Stability:	-20°C/1 year