

Antibodies Product Datasheet

Product Name:	Beta-defensin 1Antibody
Catalog Number:	C40351
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Beta-defensin 1
Concentration:	1 mg/ml
Gene Name:	DEFB1 BD1 HBD1
Protein Name:	Beta-defensin 1
Other Name:	Beta-defensin 1 (BD-1) (hBD-1) (Defensin, beta 1)
Human Swiss Prot No.:	P60022
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Beta-defensin 1
Storage and Stability:	-20°C/1 year