

Antibodies Product Datasheet

Product Name:	Beta-defensin 4AAntibody
Catalog Number:	C40352
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Beta-defensin 4A
Concentration:	1 mg/ml
Gene Name:	DEFB4A DEFB102 DEFB2 DEFB4; DEFB4B
Protein Name:	Beta-defensin 4A
Other Name:	Beta-defensin 4A (Beta-defensin 2) (BD-2) (hBD-2) (Defensin, beta 2) (Skin-antimicrobial peptide 1) (SAP1)
Human Swiss Prot No.:	O15263
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Beta-defensin 4A
Storage and Stability:	-20°C/1 year