

Antibodies Product Datasheet

Product Name:	BMP4 Antibody
Catalog Number:	C40360
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BMP4
Concentration:	1 mg/ml
Gene Name:	BMP4
Protein Name:	BMP4
Other Name:	
Human Swiss Prot No.:	Q06826.1 P21275.2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat BMP4
Storage and Stability:	-20°C/1 year