

Antibodies Product Datasheet

Product Name:	CB082 Antibody
Catalog Number:	C40403
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CB082
Concentration:	1 mg/ml
Gene Name:	C2orf82 UNQ830/PRO1757
Protein Name:	CB082
Other Name:	
Human Swiss Prot No.:	Q6UX34
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CB082 at Human/Mouse
Storage and Stability:	-20°C/1 year