

Antibodies Product Datasheet

Product Name:	CCD19 Antibody
Catalog Number:	C40419
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CCD19
Concentration:	1 mg/ml
Gene Name:	CCDC19 NESG1
Protein Name:	CCD19
Other Name:	Coiled-coil domain-containing protein 19, mitochondrial (Nasopharyngeal epithelium-specific protein 1)
Human Swiss Prot No.:	Q9UL16
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CCD19
Storage and Stability:	-20°C/1 year