

Antibodies Product Datasheet

Product Name:	CCL18 Antibody
Catalog Number:	C40429
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CCL18
Concentration:	1 mg/ml
Gene Name:	CCL18 AMAC1 DCCK1 MIP4 PARC SCYA18
Protein Name:	CCL18
Other Name:	C-C motif chemokine 18 (Alternative macrophage activation-associated CC chemokine 1) (AMAC-1) (CC chemokine PARC) (Dendritic cell chemokine 1) (DC-CK1) (Macrophage inflammatory protein 4) (MIP-4) (Pulmonary and activation-regulated chemokine) (Small-induc P55774
Human Swiss Prot No.:	P55774
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CCL18
Storage and Stability:	-20°C/1 year